



Balance TRB-605

www.labomate.com | info@labomate.com

Description

Balance TRB-605 is a compact device with high stability; the sensor uses a unique locking device with ac/dc level, an LCD backlight display, and a weighing reaction speed range of 0 to 3000 g with precision of 10 mg. It is intended to be used in conjunction with a printer or computer.

Features:

- It has ac/dc, level, LCD backlight display, fast weighing reaction time, and high stability
- The sensor employs a one-of-a-kind locking mechanism
- Counting function, percentage, automatic calibration, multiple calibration options, calibration weights
- With unit conversion, stability, memory function, animal weighing, full range gross weight reduction, and so on
- With an RS232 interface for connecting to a printer or computer

Specifications :

Accuracy	10 mg
Minimum Reading	0.01g
Pan Size	φ130 mm
Range	0 to 3000 g
Size	310×235×150 mm
Weight	1 Kg



E. Cherry St. Springfield MO 65802, USA
Email: info@labomate.com | Website: www.labomate.com