



## Constant Temperature Metal Bath TRCM-601

[www.labomate.com](http://www.labomate.com) | [info@labomate.com](mailto:info@labomate.com)



## Description

Constant Temperature Metal Bath TRCM-601 uses a simpler and more convenient temperature correction method. It can transport a wide range of modules. It has a programming number of 9 and operates at a temperature range of 5 to 40 °C. It has a patented module clamping, pressing method, and a temperature control range of 5 to 110 °C. It has an LCD screen with a touch keyboard for placing and operating parameters, and it is designed to be used in a metal bath. It can cool as well as heat.

## Features:

- The temperature adjustment range of 15 to 110 °C is excellent
- The maximum temperature accuracy is 0.5 °C
- It can transport a variety of modules
- Assist in the editing and storage of programs
- The temperature correction method is more straightforward and convenient
- The patented technique of clamping and pressing modules
- It Can be used for both cooling and heating

## Specifications :

|                                                |                                    |
|------------------------------------------------|------------------------------------|
| <b>Dimensions L×W×H</b>                        | 186×267×145 mm                     |
| <b>Display</b>                                 | LCD                                |
| <b>Input Voltage</b>                           | 100 to 120V/220 to 240V, 50Hz/60Hz |
| <b>Maximum Cooling Rate</b>                    | 5 °/min                            |
| <b>Maximum Heating Rate</b>                    | 5.5 °/min                          |
| <b>Power</b>                                   | 180 W                              |
| <b>Programming</b>                             | 9                                  |
| <b>Temperature Control Accuracy 20 To 45 °</b> | ± 0.5 °                            |
| <b>Temperature Control Range</b>               | lowest is 25 °, highest is 110 °   |
| <b>Temperature Range</b>                       | -5 ° / 110 °                       |
| <b>Temperature Uniformity 20 To 45 °</b>       | ± 0.5 °                            |
| <b>Time Setting Range</b>                      | 1min to 99 h 59 min                |
| <b>Working Environment Humidity</b>            | 0.8                                |
| <b>Working Environment Temperature</b>         | 5 to 40 °                          |



E. Cherry St. Springfield MO 65802, USA  
Email: [info@labomate.com](mailto:info@labomate.com) | Website: [www.labomate.com](http://www.labomate.com)