



## Mini Vertical Electrophoresis cell TRMV-601

[www.labomate.com](http://www.labomate.com) | [info@labomate.com](mailto:info@labomate.com)

## Description

Mini Vertical Electrophoresis cell TRMV-601 has a capacity of 1-2 pieces and it is used for pre-casting and hand-casting gel. There are no more than four pieces of gel available at any given time. It has a multifunctional gel scraper built-in. The superior comb does not slow gel polymerization. The gel is separated from the air by a built-in kit. Ensures a high level of leak resistance. Supplied with an 8373mm gel size.

## Features:

- The filler strip is stable to the glass to avoid leakage of gel
- A built-in kit separates the gel from the air
- The gel-making kit safeguards accurate aligning
- The distinctive comb prepares not to stop gel polymerization
- Excellent leak proof ness
- The thickness and well-marked comb are understandable
- Equipped with a multi-functional gel scraper

## Specifications :

<b>Comb Specification</b>	1.0mm: 10 and 15 wells
	0.75mm: 10 and 15wells
	1.5mm: 10 and 15 wells
<b>Dimension (L ×W×H)</b>	180×130×160mm
<b>Gel Dimension</b>	83×73mm
<b>Quantity</b>	1 to 2 piece(s)
<b>Short glass Plate</b>	100×73mm
<b>Spacer glass plate</b>	100×83mm
<b>Thickness</b>	1.0mm: 10, 15 wells (standard), 0.75mm: 10 & 15 wells, 1.5mm: 10 & 15 wells (optional)



E. Cherry St. Springfield MO 65802, USA  
Email: [info@labomate.com](mailto:info@labomate.com) | Website: [www.labomate.com](http://www.labomate.com)