



Ultrasonic Cleaner TRUC-605

www.labomate.com | info@labomate.com

Description

Ultrasonic Cleaner TRUC-605 has an ultrasonic vibration frequency of interaction composed of high-quality SUS304 stainless steel stampings. It has a 15 L capacity, a heating temperature range of 20 to 80°C, and a digital timer with a setting range of 0 to 99 minutes. It does not have a drainage ball valve.

Features:

- The intensity of ultrasonic vibration, thoroughly solve the ultrasonic cleaning
- Its work-performed tank is made of high-quality SUS304 stainless steel molding
- The bladder is made of high-quality steel material, and the cladding is made of aluminum
- Electronic timer with a period range of 0 to 99 minutes
- The drainage ball valves are installed in drainage systems 6L and higher models

Specifications :

1 Unit Packing	46×41×43 cm
1 Unit Weight	10 Kg
Dimensions	36×33×31 Cm
Drainage	With
Heating Power	400W
Heating Temperature	20 To 80?
Overall Packing	88×49×79 Cm
Rinse Tank Size	33×30×15 cm
Timer	0 To 99min
Ultrasonic Frequency	28/40khz
Ultrasonic Power	360W
Volume	15 L
Weight	42 Kg



E. Cherry St. Springfield MO 65802, USA
Email: info@labomate.com | Website: www.labomate.com